

BIOMEDICA & ENVIRONMENTAL MONITORING



Products

AirSampler

- Portable Air Samplers
- Compressed Gases
- Isolator Solutions

Particle Counter

- Stand-alone Counters
- External Vacuum Pump
- Plug&Play Systems

Environmental Monitoring Software

- Truly Paperless
- Full Traceability
- User Friendly



PRODUCT INFORMATION

AirSampler



Microbiological control of air and gases

- Standard (1 – 3 heads)
- Space-saving: suitable for isolators/RABS
- TrioGas for compressed gases
- FDA 21 CFR part 11 compliant, USP/EP/JP conform
- IQ/OQ/PQ
- Full data integrity (BAS Software)
- Daily shift heads (disposable)
- Cable or induction charging
- Remote aspirating systems
- Calibration service
- Wide range of accessories

Particle Counters



Stand-alone, mobile or remote-control instruments

- Plug & Play instruments
- 21 CFR part 11 compliant, ISO 14644 conform, EUGMP, Annex 1
- Stainless steel housing
- H2O2 resistant
- 8 configurable particle size channels (0,2-10µm)
- Data export via USB, csv-file or network
- Communication via NIOBus- or Modbus-protocol

Environmental Monitoring Software



Software for Reliable and Effective Environmental Monitoring

- User-friendly, truly paperless
- FDA 21 CFR Part 11 and EU GMP Annex 11 compliant
- Cloud-based – secure and already validated
- Full data integrity compliance
- Quick implementation
- Audit trail
- Planning and data analysis developed for users

