

## BIOMEDICA & READY PREPARED MEDIA



### Our Partner

**bioTRADING Benelux B.V.** manufactures ready prepared culture media since 1979 with over 40 years of experience and is accredited by the Dutch Council of Accreditation (reg.nr. L588) in accordance with the recognized International Standard ISO/IEC 17025 'General Requirements for the competence of testing and calibration laboratories'. The accredited methodology conforms to ISO 11133. It is the policy of bioTRADING to achieve the highest degree of product quality, ensuring that processes and services are continuously improved and are also in line with international standards, FDA/ cGMP/ GLP guidelines and customer requirements. The ultramodern production facility in Mijdrecht, The Netherlands, that opened in 2009 is designed and equipped to meet pharmaceutical regulations, for Ready Prepared Culture Media.

The facility and environmental monitoring (including microbiological and particle measurement) systems are 'FDA 21 CFR part 11 compliant' and manufacturing takes place in cleanrooms class 5 to 8 (conform ISO14644-1).

Today, bioTRADING offers more than 1.000 different Ready Prepared Culture media recipes. Custom specialised media can be manufactured to order in both small or large batch volumes, to a preferred recipe, bottle or tube size, and volume. Current batches range in size from small quantities of 10 units, to batch sizes of over 25.000.



## PRODUCT INFORMATION

### Standard Pre-Poured Media



- Petri Dish available in 55mm, 90mm, 100mm and 140 mm
- Square Petri Dish available in 120mm
- Contact Plates available in 30mm and 55mm
- Single, double & triple wrapped
- Extended shelf life, deep fill & ambient storage options
- Over 1.000 formulations

### Bottled Media



- Available in Alpha, Infusion & Polypropylene Bottles
- From 30 to 1000 ml
- Complete range of bottled agars & broths
- Options for screw, septum und felseap lids

### Infusion Bags & Vial



- Infusion bags from 100 ml to 5000 ml
- Infusion bags available with septum, luer lock connection or CPC
- Vials with septum or open felseap
- From 10 to 100 ml – all the same passage size
- Available empty or filled with TSB
- Dehydrated and concentrated TSB available

### Tubes



- Snap cap tubes available from 15 to 65 ml
- Screw cap tubes from 15 to 23 ml
- Plastic tubes from 5 to 20 ml

### QBigSet



The media bottles are packed in a plastic box, wrapped in a Tyvek-header-bag and then sterilised using Eto gas sterilisation. The unique solution results in bottles that are free from contamination and can be transferred immediately into your cleanroom or isolator. It is possible to configure the kit to your exact requirements, be it one type of bottled media per box or a mix. You can also specify the number of wraps for each box.

### Custom Solutions

- Manufactured to order in preferred batch sizes
- Any media, any format

